de	Gebrauchsanweisung 📄 3	no	Bruksanvisning 📄 60
en	Instructions for Use 📄 10	pl	Instrukcja obsługi 📄 67
fr	Notice d'utilisation 📄 17	ru	Руководство по эксплуатации 📄 74
es	Instrucciones de uso 📄 25	sr	Uputstvo za korišćenje 📄 81
pt	Instruções de uso 📄 32	sl	Navodila za uporabo 📄 88
it	Istruzioni per l'uso 📄 39	sk	Návod na použitie 📄 95
da	Bruksanvisning 📄 46	cs	Návod k použití 📄 102
fi	Käyttöohje 📄 53	tr	Kullanım talimatları 📄 109

Sauerstoff-Druckgasflaschen Compressed Oxygen Cylinders

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1 For your safety

1.1 General safety statements

Before using this product, carefully read these Instructions for Use and those of the associated products.

Strictly follow the Instructions for Use. The user must fully understand and strictly observe the instructions. Use the product only for the purposes specified in the Intended use section of this document.

Do not dispose of the Instructions for Use. Ensure that they are retained and appropriately used by the product user.

Only trained and competent users are permitted to use this product. The operator is responsible for ensuring that the product is used correctly.

Naked flames, light and smoking are forbidden whilst working with compressed air cylinders.

Follow the local and national guidelines pertaining to this product. Observe the regulations concerning operation and maintenance of compressed air cylinders.

Only trained and competent personnel are permitted to inspect, repair and service the product. Dräger recommend a Dräger service contract for all maintenance activities and that all repairs are carried out by Dräger.

Use only genuine Dräger spare parts and accessories, or the proper functioning of the product may be impaired.

Do not use a faulty or incomplete product. Do not modify the product. Faults and damage to compressed air cylinders or safety mechanisms must be reported to the operator without delay.

Notify Dräger in the event of any product or component fault or failure.

1.2 Definitions of alert icons

The following warning symbols are used in this document to provide and highlight areas of the associated text that require a greater level of awareness from the user. The meanings of the symbols are as follows:

WARNING

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, could result in physical injury, or damage to the product or environment. It may also be used to warn against unsafe practices.

NOTICE

Indicates additional information on how to use the product.

2 Description

2.1 Product overview

The body of the compressed oxygen cylinder is made of carbon fibre reinforced plastic, steel or aluminium. The compressed oxygen cylinder is equipped with a cylinder valve. The cylinder valve is available with different connections.

The period of service depends on the volume of the compressed oxygen cylinder and the breathing rate of the carrier.

2.2 Intended use

The compressed oxygen cylinder is used in approved closed-circuit breathing apparatus (e. g. Dräger PSS BG 4). It provides the user with an independent supply of breathable air in a harmful, contaminated or oxygen-deficient atmosphere.

2.3 Approvals


Compressed oxygen cylinders with a valve are pressure vessels. They were developed for use in closed-circuit breathing apparatus. They meet the requirements specified in the Pressure Equipment Directive. They have a design approval and were either tested individually or through production monitoring by the notified body, TÜV Nord (technical inspection authority).

NOTICE

Compressed oxygen cylinders must always be used with a valve that is approved for the cylinder otherwise the approval will expire.

2.4 Explanation of type-identifying marking and symbols

00721703.eps

- 1 ID No. of notified body
- 2 Symbol "Follow instructions for use"
- 3 Hazard symbol: "Hazardous material class 2.2 - non-flammable/combustible gases"
- 4 Hazard symbol: "ADR class 5.1 - materials with a combustible (oxidising) effect"
- 5 Warning: "Oxygen compressed, send back with residual pressure"
- 6 Pressure and temperature specifications, intended use: A
- 7 Additional information

The illustrated name plate is an example. The specifications may vary depending on the compressed oxygen cylinder.

Compressed oxygen cylinders made of carbon fibre reinforced plastic have a label with the type identifying marking.

Compressed oxygen cylinders made of steel or aluminium have the serial number, test date and reference number of the notified body stamped on the shoulder. Additional tests of these compressed oxygen cylinders are documented on labels.

"Intended use: A" means that the compressed oxygen cylinder can be used in closed-circuit breathing apparatus.

"PS": Operating pressure

"PT": Test pressure

"TS": Operating temperature

3 Use

3.1 Information on use

- Do not completely empty the compressed oxygen cylinder as moisture could get in and cause corrosion or contamination inside. Dräger recommends leaving a minimum pressure of 2 bar in the compressed oxygen cylinder. Compressed oxygen cylinders that have been completely emptied by mistake must be dried and tested. In this case, send the compressed oxygen cylinder to Dräger to be checked.
- Do not hold the thumb wheel on the valve to carry the compressed oxygen cylinder. Hold it by the valve housing or the cylinder body.
- Compressed oxygen cylinders can be heavy and thus difficult to handle. Improper handling of compressed oxygen cylinders can lead to sprains, over-exertion, falls, bruising or bone fractures. Avoid accidents by taking sufficient precautionary measures.
- Only open or close the valve on the compressed oxygen cylinder by hand. Do not use tools! Do not use force to turn the valve on the compressed oxygen cylinder.
- The compressed oxygen cylinder contains pressurised gas. It can explode if it is heated.
- Escaping gas can cause a fire or intensify it. Keep the compressed oxygen cylinder away from flammable materials!

3.2 Preparation for use

3.2.1 Filling the compressed oxygen cylinder

WARNING

Explosion hazard!

Compressed, flowing oxygen has an oxidising effect and causes oil and grease to ignite immediately.

The inside of the cylinder, the valve thread and the sealing elements must be free from oil and grease.

Compressed oxygen cylinders may only be filled when they:

- have a label or stamp with the test date and identifier of the accredited inspection authority and the specified test interval has not been exceeded,
- do not show any signs of damage that could be dangerous (e. g. a damaged valve housing, thumb wheel or a leaky valve),
- are in good condition and show no signs of damage, corrosion or moisture on the valve,
- are in accordance with national guidelines.

Compressed oxygen cylinders may only be filled with oxygen of the following quality characteristics:

- Purity: ≥ 99.5 vol. % O_2
- Water content: ≤ 0.05 mg/L oxygen
- Tasteless and odourless
- The maximum permissible contaminants are subject to national regulations, in Germany, for example:
Inert gases ($N_2 + Ar$): 5000 ppm
Carbon monoxide (CO): 5 ppm
Carbon dioxide (CO_2): 300 ppm
Oil: ≤ 2.5 mg/m²

The compressed oxygen cylinders must be filled at the correct nominal filling pressure as specified on the cylinder body and valve housing.

WARNING

Danger of bursting!

Do not overfill the compressed oxygen cylinders! They do not have a safety valve and could burst.

To avoid overfilling, Dräger recommends the use of a pressure relief valve on the oxygen filling facility.


The temperature of the oxygen and compressed oxygen cylinder can increase during the filling process. This can lead to an incomplete filling of the cylinder. To avoid this effect to a large extent, Dräger recommends a filling pressure increase of 27 bar/minute. Check the cylinder pressure at ambient temperature after the filling process and top up the compressed oxygen cylinder, if required.

Observe the instructions of the filling station manufacturer.

Compressed oxygen cylinders made of carbon fibre reinforced plastic can develop cracks in the paint coating at the cylinder neck or base when the cylinders are being filled. This is due to the varying distention of the composite material on the shell and at the cylinder neck/base. This does not impair the operational safety of the compressed oxygen cylinders.

3.2.2 Installing the compressed oxygen cylinder

1. Make sure the cylinder valve is closed.
2. Remove the lock nut from the valve outlet.
3. Place the compressed oxygen cylinder at the designated position in the closed-circuit breathing apparatus; secure and connect it¹⁾.
4. Carry out a high-pressure leak test in accordance with the Instructions for Use for the closed-circuit breathing apparatus.

3.3 During use

The compressed oxygen cylinder is used in accordance with the Instructions for Use for the closed-circuit breathing apparatus.

3.4 After use

3.4.1 Removing the compressed oxygen cylinder

1. Make sure the cylinder valve is closed.
2. Remove the compressed oxygen cylinder from the closed-circuit breathing apparatus²⁾.
3. Screw the lock nut onto the valve outlet.

CAUTION

Risk of corrosion!

Make sure the valve is closed. Otherwise moisture can get into the compressed oxygen cylinder and cause contamination or corrosion.

3.4.2 Checking the compressed oxygen cylinder

- If the compressed oxygen cylinder has come into contact with corrosive media (e. g. acid), it has to be depressurised and Dräger must be consulted to clarify the subsequent course of action.
- If the compressed oxygen cylinder had to withstand a heavy impact because it was, e.g. dropped, then it must be depressurised and sent to Dräger for assessment.
- Compressed oxygen cylinders must be discarded in the following cases:
 - The compressed oxygen cylinder has permanent deformations such as dents in the material, grooves or pulping that are longer than 20 mm or deeper than 0.25 mm.
 - The compressed oxygen cylinder shows signs of the paintwork or composite material flaking or blistering or other structural damage.
 - The compressed oxygen cylinder shows signs of abrasive wear.

¹⁾ Refer to the Instructions for Use for the closed-circuit breathing apparatus

²⁾ Refer to the Instructions for Use for the closed-circuit breathing apparatus

- The cylinder material, label or paint are burnt or blackened or the paint has melted.
- A clear allocation of the compressed oxygen cylinder is not possible because the label is illegible.
- If a clear-cut assessment of the damage is not possible, then Dräger should be consulted to clarify the subsequent course of action.
- If the date for the periodic test has been exceeded, the compressed oxygen cylinder may not be filled or used. A periodic test must be arranged.

4 Maintenance

4.1 Maintenance intervals

Work to do	Before use	After use	Monthly	Every 5 years
Check filling pressure	X		X	
Fill the compressed oxygen cylinder ¹⁾		X		
Periodic test ²⁾				X ³⁾
Major overhaul of the valve ⁴⁾				X


- 1) see section 3.2.1 on page 13
- 2) by an accredited inspection authority
- 3) observe national guidelines
- 4) by Dräger or a service technician trained by Dräger

5 Transportation

Secure compressed oxygen cylinders against falling over and slipping, and against heat, chemicals and fire.

For transportation, leave a residual pressure of 2 bar in the compressed oxygen cylinder.

Make sure that the lock nut closes the valve connection.

6 Storage

The valve must remain closed during storage. The valve must not remain open when the compressed oxygen cylinder is fully depressurised. Dräger recommends leaving a minimum pressure of 2 bar in the compressed oxygen cylinder. Dirt, dust and moisture get into empty compressed oxygen cylinders with an open valve. In this case, send the compressed oxygen cylinder to Dräger to be checked.

Only store compressed oxygen cylinders in designated, well ventilated storage areas. Full and empty compressed oxygen cylinders must be stored separately.

Check stored compressed oxygen cylinders regularly (at least once a month) for signs of leakage or other problems.

7 Disposal

Compressed oxygen cylinders made of steel or aluminium can be used without a time limit.

The service life of compressed oxygen cylinders made of carbon fibre reinforced plastic can be limited. Compressed oxygen cylinders made of carbon fibre reinforced plastic with a service life limit are identified by FIN or FINAL and the date of expiry. Do not use these compressed oxygen cylinders after the service life has expired!

Compressed oxygen cylinders made of carbon fibre reinforced plastic without a service life limit are marked NLL.

NOTICE

The stipulated service life is only attained if the compressed oxygen cylinders are used correctly and checked regularly.

You can return this product to Dräger at a charge. A decision on its subsequent use will be made when the product is returned. Please contact your national Dräger Sales Organisation or Dräger for more information.

Please send compressed oxygen cylinders in approved packaging to the following address:

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8 Technical data

Operating pressure	200 bar
Test pressure	300 bar
Operating temperature	-30 °C ... +60 °C
Storage temperature (for filled compressed oxygen cylinders)	-30 °C ... +60 °C
Connection for cylinder valve	W21,8x1/14" M24x2 (according to EN 144) G3/4" (according to ISO 228) W22x1/14"

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© Dräger Safety AG & Co. KGaA
Dräger Safety AG & Co. KGaA
TÜV NORD CERT GmbH